

AuraRing Pro Gen 3

Customer Support Knowledge Base

1. Product Overview & Specifications

1.1 Product Identity

- **Product Name:** AuraRing Pro Gen 3
- **Category:** Smart Health & Wellness Wearable (Smart Ring)
- **Model Numbers:** AR3-TIT-BLK (Matte Black), AR3-TIT-SLV (Brushed Silver), AR3-TIT-GLD (Rose Gold) • **Target Audience:** Fitness enthusiasts, sleep & stress tracking professionals, and users seeking discreet biometric monitoring without screen distractions.

1.2 Technical Specifications Table

Feature / Component	Specification
Material & Build	Aerospace-grade Titanium with hypoallergenic medical-grade PEEK inner molding
Weight & Dimensions	4 to 6 grams (depending on size); Width: 7.9 mm; Thickness: 2.5 mm
Sizing	Standard US Sizes 6 through 13 (no half sizes)
Sensors	Dual PPG Optical Pulse Sensors (Red, Green, IR), Skin Temperature Sensor (PTC Thermistor, ±0.1°C), 3D Accelerometer, Optical SpO2 Sensor
Battery Life	Up to 7 days on a single charge (varies based on continuous SpO2 frequency)
Battery Capacity	18–22 mAh (varies by ring size) Lithium-Polymer
Charging Method	Magnetic Contact Induction Cradle (USB-C input, 5V/1A)
Charge Time	20% to 100% in 55 minutes; 0% to 100% in 75 minutes
Water Resistance	10 ATM (up to 100 meters / 330 ft); safe for swimming, showers, sauna up to 60°C
Connectivity	Bluetooth Low Energy 5.3 (BLE) with AES-128 hardware encryption
Onboard Memory	Stores up to 14 days of biometric data offline without smartphone syncing
Compatibility	iOS 16.0+; Android 11.0+; Apple Health & Google Health Connect integration

2. Setup & Getting Started

2.1 Sizing and Fit Guidelines

- **Sensor Alignment:** The small raised bumps on the inner band house the optical PPG sensors. These sensors must sit on the palm side (underside) of your finger for accurate readings.
- **Finger Recommendation:** The index finger of either hand provides the highest physiological signal accuracy. Middle and ring fingers are acceptable secondary alternatives.
- **Fit Check:** The ring should feel snug without shifting during active hand motion, but must slide smoothly over your knuckle without forcing or soapy water.

2.2 First-Time Pairing Process

1. Connect the **AuraRing Pro Charging Dock** to a certified 5V/1A or 5V/2A USB-C power source. AuraRing Health Inc. —
2. Place the AuraRing Pro into the dock. The dock LED will turn solid **Amber**, confirming active charging. 3. Download the **Aura Health Companion App** from the Apple App Store or Google Play Store. 4. Enable Bluetooth on your smartphone and keep the device within 1 meter of the dock.
5. Open the app, select **Set Up New Device**, and follow on-screen instructions.
6. Confirm the 4-digit pairing code shown in the app matches the label underneath your dock.
7. Keep the ring docked until any initial mandatory over-the-air firmware updates conclude.

3. Battery, Charging & LED Indicators

3.1 Charging Dock LED Reference Table

LED Status	Meaning	Required Action
Pulsing Amber	Ring is docked and currently charging (< 90%).	Leave ring on cradle until charge completes.
Solid Green	Ring is docked and fully charged (100%).	Ring is ready for wear.
Blinking Red (3 Hz)	Docking misalignment or foreign debris on contacts.	Remove ring, clean gold contact pins with dry microfiber cloth, and reseal.
Pulsing Blue	Dock is in Bluetooth pairing / discovery mode.	Complete pairing flow inside the Aura Health App.
Solid White	Dock firmware update in progress.	Do not disconnect power cable until LED shuts off.

3.2 Battery Optimization & Care

- Avoid leaving the ring battery fully depleted (0%) for more than 14 consecutive days to preserve long-term chemical cell integrity.
- Optimal charging temperature is 15°C to 30°C (59°F to 86°F). Internal thermal cutoffs automatically suspend charging if temperature exceeds 45°C (113°F).

4. Troubleshooting & Error Codes

4.1 System Error Codes Table

Error Code	Description	Root Cause	Resolution Steps
ERR-SYNC-101	Data Sync Timeout	Bluetooth transfer interrupted during packet transmission.	1. Toggle phone Bluetooth. 2. Force close and reopen app. 3. Place ring on powered dock to trigger wake-up packet.
ERR-SENS-204	Biometric Signal Degradation	Optical sensors cannot register pulsatile capillary flow.	1. Verify sensor bumps face palm side. 2. Clean inner window with 70% isopropyl alcohol. 3. Ensure proper ring size.
ERR-FW-302	Firmware Checksum Error	OTA firmware packet corrupted during wireless transfer.	1. Keep ring docked with >50% charge. 2. Clear app cache. 3. Restart OTA update from Settings > Device.
ERR-BATT-401	Rapid Drain Protection	Background continuous sampling sensor loop hang.	Perform a Soft Hardware Reset (Section 4.3).

4.2 Sync & Connectivity Troubleshooting

If the app reports '**Ring Disconnected**':

1. Verify battery level by docking the ring. If dock LED does not light up, charge for at least 20 minutes before proceeding.
2. Confirm Bluetooth and Nearby Device permissions are granted in phone OS settings.
3. **For Android:** Confirm Location Services are enabled (required for Android BLE beacon discovery).
4. Remove 'AuraRing Pro' from your phone's bonded Bluetooth devices, then re-pair through the app.

4.3 Reset Procedures

- **Soft Hardware Reset (Retains Cached Data):** Seat ring in powered dock. Use a SIM pin to press and hold the pinhole button under the dock for **5 seconds** until the dock LED flashes green twice.
- **Factory Reset (Erases Local Memory):** Open app > Profile > Settings > Device Settings > Unpair & Factory Reset. Then hold dock pinhole button for **15 seconds** until LED pulses Red 3 times.

5. Maintenance, Care & Environmental Limits

5.1 Cleaning Protocols

- Wash ring weekly using mild liquid hand soap and lukewarm water.
- Dry thoroughly with a clean microfiber cloth before placing on the charging cradle.
- **Prohibited:** Ultrasonic jewelry baths, acetone, pure chlorine bleach, and abrasive metal polishes.

5.2 Water & Thermal Limits

- **Approved:** Swimming pools, ocean salt water (rinse with freshwater post-swim), showers, handwashing.
- **Temperature Range:** Operating range is -10°C to 50°C (14°F to 122°F). Exposure to dry saunas >70°C or ice baths <-15°C activates protective automatic thermal shutdown.

6. Warranty, Returns & Support Policies

6.1 Return Window & Sizing Exchanges

- **30-Day Money-Back Trial:** Returnable within 30 calendar days of confirmed delivery in original packaging for a full refund.
- **Free Size Exchange:** One complimentary sizing exchange is included within the first 30 days of receipt.

6.2 Warranty Coverage

- **Duration:** 1-Year Limited Hardware Warranty (2 years for European Union and UK residents).
- **Covered:** Battery capacity loss below 50% under normal usage within warranty term; optical sensor malfunction; BLE chipset failure; separation of titanium shell from inner resin molding.
- **Excluded:** Normal cosmetic wear (scuffs, knurling marks from gym weights); submersion beyond 100m depth; damage caused by non-certified third-party charging docks.

6.3 Support Contact Channels

- **In-App 24/7 Live Chat:** Aura Health App > Profile > Help & Support
- **Email Support:** support@auraringhealth.example.com (Average turnaround < 4 business hours)
- **Help Desk Portal:** <https://help.auraringhealth.example.com>

7. Frequently Asked Questions (FAQ)

Q: Can I wear my AuraRing Pro while lifting weights or using barbells?

A: While titanium is durable, knurled steel barbells, dumbbells, and kettlebells can cause cosmetic scratching and ring pinch. We advise wearing a protective silicone ring guard or removing the ring during heavy barbell workouts.

Q: Does the AuraRing Pro require an ongoing subscription to view health data?

A: No. Core tracking metrics (Sleep Score, Resting Heart Rate, HRV, SpO2, Daytime Activity) are free forever. An optional Aura Health+ subscription (\$5.99/mo) offers multi-year longitudinal biomarker analysis and AI sleep coaching.

Q: Can I pair the ring to multiple smartphones simultaneously?

A: No. For security and biometric data privacy, the ring maintains a single active BLE binding. To migrate phones, sync existing data to the cloud, unpair in settings, and pair on the new phone.

Q: What happens if I travel without my phone?

A: The onboard non-volatile memory logs up to 14 days of full biometric records. When you reconnect with your phone, all logged telemetry synchronizes automatically.